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**JOYTECH Healthcare Co., Ltd.**

Adresse : No.365, Wuzhou Road, 311100 Hangzhou,  
Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA

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## USER'S MANUAL

### Warnings

- ⚠ Read instructions thoroughly before using digital thermometer.
- ⚠ Choking Hazard: Thermometer cap and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision.
- ⚠ Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- ⚠ Do not place thermometer battery near extreme heat as it may explode.
- ⚠ Remove battery from the device when not in operation for a long time.
- ⚠ The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
- ⚠ Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- ⚠ Do not bite the thermometer. Doing so may lead to breakage and/or injury.
- ⚠ Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- ⚠ After each use, disinfect the thermometer especially in case the device is used by more than one person.
- ⚠ Do not force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
- ⚠ Do not use thermometer orally after being used rectally.
- ⚠ For children who are two years old or younger, please do not use the devices orally.
- ⚠ If the unit has been stored at temperatures over 5°C - 40°C, leave it in 5°C - 40°C ambient temperature for about 15 minutes before using it.
- ⚠ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- ⚠ PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the [ME EQUIPMENT or ME SYSTEM], including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.
- ⚠ It is not intended for use in the oxygen rich environment and presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- ⚠ Do not put the thermometer in direct sunlight or with cotton wool, otherwise the accuracy will be affected.
- ⚠ The equipment should not be cleaned and disinfected while in use.
- ⚠ Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/ or patient is established.

### INTENDED USE

The digital thermometers are intended to measure the human body temperature in regular mode orally, rectally and under the arm. And the devices are reusable for clinical or home use on people of all ages, including children under 8 years old with adult supervision.

### INTENDED USER/INDICATION FOR USE

Patient or healthy people who want to measure their body temperature, and medical workers who want to measure the patient's body temperature.

Patient target group :Patients or healthy people who need or want to have their temperature taken.

### PLEASE READ CAREFULLY BEFORE USING

The predictive digital thermometer provides a quick and highly accurate reading of an individual's body temperature. Predictive-read thermometers are quicker than actual read thermometer. Predictive read type thermometers display the temperature results in a short period of time which are equivalent to balance temperature after 5 minutes according to the particular algorithm. Therefore, users only need about 5 seconds to take the temperature readings. Due to measurement sites of predict thermometer are different, the reading time also may be different. (See below Figure 1).


Figure 1

### CONTENTS

1 thermometer - 1 user's manual- 1 Storage Case

### PRODUCT ILLUSTRATION


Figure 2

### PRÉCAUTIONS

→ The performance of the device may be degraded should one or more of the following occur:

- Operation outside the manufacturer's stated temperature and humidity range.
- Storage outside the manufacturer's stated temperature and humidity range.
- Mechanical shock (for example, drop test) or degraded sensor.
- Patient temperature is below ambient temperature.

→ Portable and mobile RF communications can affect the device. The device needs special precautions regarding EMC according to the EMC information provided in the accompany documents.

→ Do not use the devices in the MR environment.

### SYMBOL EXPLANATION

Symbol	Description	Symbol	Description
	Caution ! Read the provided documents (IFU)		Date of manufacture
	Refer to instruction manual/ booklet		Storage and Transportation Temperature Limit
	Type BF applied part		Storage and Transportation Humidity limitation
	Batch Code		Atmospheric pressure limitation
	Authorized representative in the European Community		Medical device
	Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.		Reusable packaging and notice
	CE conformity marking 0123 represent Notified Body identification number		MR Unsafe
	Manufacturer		The first num.2:Protected against solid foreign objects of Ø 12.5 mm and greater. The second num.7:Protected against the effects of temporary immersion in water.
	Imported by		Distributed by

### SPECIFICATIONS

Type	Digital Thermometer (Predictive)			
Measure Range	28,0°C - 42,9°C			
Accuracy	±0.1°C during 35.5°C-42.0°C at 18°C-28°C ambient operating range ±0.2°C for other measuring and ambient operating ranges			
Operating mode	Adjusted mode: Oral mode/Rectal mode/Underarm mode Direct mode: Bath mode			
Display	Liquid crystal display, 3 1/2 digits			
Memory	The last ten memories			
Battery	One 3.0V DC button battery type CR2032			
Battery life	Approx. 200 hours of continuous operation or 1 year with 3 measurements per day			
Dimension	14,3 cm × 2,5 cm × 1,4 cm (L x W x H)			
Weight	Approx. 20 grams including battery			
Expected service life	Three years			
Ambient operating range:	Temperature : 5°C - 40°C - Hygrometry : 15 % - 95 % HR Atmospheric Pressure : 70 kPa - 106 kPa			
Storage and transportation condition	Temperature : -20°C - 55°C - Hygrometry : 15 % - 95 % HR Atmospheric Pressure : 70 kPa - 106 kPa			
Clinical precision		Oral Use	Rectal Use	Armpit Use
	Clinical bias	0,09	0,06	0,07
	Agreement limits	0,72	0,71	0,71
	Clinical repeatability	0,11	0,12	0,12
Ingress Protection Rating:	IP 27			
Classification	Type BF			
Contraindication	No contraindication			

## MEMORY MODE

1. Power on the thermometer.
2. During the last memory displays, press and hold the On/Off button until entering into Memory Mode.
3. Press the button again to cycle check the last 10 memories.
4. Hold down the on/Off button for 3 seconds or wait for one minute to exit the memory mode.

## SELECT MODE

1. Power on the thermometer.
2. During the mode displays, press and hold the On/Off button until mode flashing.
3. Press the button again to cycle through the four mode options (see figure 3).
4. When the desired mode displays on screen, wait for auto entering into measurement mode.

Figure 3

## DIRECTIONS

1. Press the On/Off Button next to LCD display. A tone will sound as the screen shows (See Figure 4), followed by last recorded temperature. After showing the mode of measurement (See Figure 3), the thermometer will enter into the testing mode (see Figure 5).

Figure 4

Figure 5

2. Position thermometer in desired location (mouth, rectum, or armpit.)
  - a) **Oral Use:** Place thermometer under tongue as indicated by position « √ » shown in Figure 6. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.
  - b) **Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2) into rectum.
  - c) **Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required.
  - d) **Bath Mode:** It will display --- and then °C is still flashing (See Figure 5), then place probe tip into bath water.


Figure 6

3. The three dashes (---) will flash sequentially throughout the testing process, and at the same time the « °C » are flashing (See Figure 7). When flashing stops an alarm will beep for approximately 5 seconds. The predictive measured reading will appear on the LCD simultaneously (For example 37.0°C - see Figure 8). The measurement time of the oral/underarm/rectal mode varies with the individual.

**Note:** The minimum measurement time of Bath Mode until the signaling tone (beep) must be maintained without exception.


Figure 7

Figure 8

Figure 9

4. The thermometer will enter into the **Actual Measurement** after 3 minutes **except Bath Mode**. At the same time you can hear ten beeps as a cue (The « °C » will disappear and the Actual measurement temperature display see Figure 9), users can take the true body temperature at this time. Users can compare the predictive temperature result with true temperature result. In order to achieve better body temperature measurement result, recommend to keep the probe in mouth and rectum about 2 minutes, or in armpit about 5 minutes regardless of the beep sound and at least 30 seconds measurement intershould be maintained.
5. To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.
6. Store the thermometer in its protective case.

**Note:** Normally the buzzes are "Bi-Bi-Bi-"; Alarm beeps more rapidly when temperature reaches 37.8°C or higher, and the buzzes are "Bi-Bi-Bi-----Bi-Bi-Bi-----Bi-Bi-Bi".

## DÉPANNAGE

Error message	Problem	Solution
- - -	Temperature taken is lower than 32.0°C	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
H <sub>i</sub>	Temperature taken is higher than 42.9°C	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
Er 1	The thermometer is moved when measuring.	Do not move, and then repeat measurement.
Er 2	The predictive measurement time is greater than 15 seconds.	The ambient temperature may be too low, please measure at a higher temperature. The person may have just had a cold drink, repeat the measurement and wait 15 minutes. The thermometer is moved during the measurement. Hold still and repeat the measurement.
Other error messages	The thermometer is not working properly.	Remove the battery, wait 1 minute and reinsert it. If the message reappears, contact your dealer for repair.
🔋	Dead battery : Only the battery icon appears, no measurement can be taken.	Replace battery.

## BATTERY REPLACEMENT

1. Replace battery when the symbol « 🔋 » appears in the lower right corner of the LCD display.
2. Pull battery cover off as shown in Figure 10.
3. Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2") (see Figure 11)
4. Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 3.0V DC button type CR2032. Be sure battery is installed with « - » polarity facing up. (see Figure 12)
5. Slide battery chamber back into place and attach cover.

Figure 10

Figure 11

Figure 12

## CLEANING AND DISINFECTION

1. Immerse the thermometer probe in distilled water for at least 1 minute;
2. Using a clean, soft cloth to wipe the thermometer down to remove any residue;
3. Repeat step 1 and 2 for three times until no soil is seen with visual inspection after cleaning;
4. For thoroughly clean and disinfection, please use method A or B:
  - Method A (High level disinfection):** immerse the thermometer probe in 0.55% OPA(O-Phthaldehyde), such as CIDEX OPA, for at least 12 minutes under temperature at 20°C;
  - Method B (Low level disinfection):** Using a clean soft cloth dipped in 70% medical alcohol, wipe the probe 3 times, at least one minute for each time.
5. Repeat step 1 to 3 to remove OPA residuals;
  - Note 1:** Rectal use is not recommended for home use as OPA will not be readily available outside of a hospital. If rectal measurement is necessary, we strongly recommend high level disinfection.
  - Note 2:** Please operate according to the manual of OPA for reference.

To prevent damage to the thermometer please note and observe the following:

- Do not use benzene, paint thinner, gasoline or other strong solvents to clean the thermometer.
- Do not attempt to disinfect the sensing probe (tip) of the thermometer by immersing in alcohol, OPA or in hot water (water over 50°C) for long time.
- Do not use ultrasonic washing to clean the thermometer.

## CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used periodical, readjustment is not necessary. However, we recommend to check the calibration every two years, or whenever the thermometer's clinical accuracy is called into question. Set the temperature of the high-precision water bath to 33.00°C, 37.00°C, 42.00°C. Switch on the thermometer and insert it into the high-precision water bath, then check the thermometer's laboratory accuracy. Repeat the previous measurement process after a 10-minute stability period. Please send the complete unit to the dealer or manufacturer. The above recommendations do not replace legal requirements. The user must always comply with the legal requirements regarding measurement control, functionality and accuracy of the device as required under the relevant laws, directives or ordinances where the device is used.

## LIMITED WARRANTY

The thermometer is guaranteed for two years from the date of purchase. If the thermometer does not function properly due to defective components or poor workmanship, we will repair or replace it free of charge. All components are covered by this warranty excluding the battery. The warranty does not cover damages to your thermometer due to improper handling. To obtain warranty service, an original or copy of the sales receipt from the original retailer is required.

This appliance conforms to the following standards:

ISO 80601-2-56 Medical electrical equipment —Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement,  
EN 60601-1-11 Medical electrical equipment —Part 1-11: General requirements for basic safety and essential performance

Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1 (Safety) standards. And the manufacturer is ISO 13485 certified.

Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.

The product is in compliance with the requirements of MDR(EU)2017/745, «0123» is the identification number of notify body.

## ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard CEI 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regards to EMC which must be published in the IFU. Portable and mobile communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

**Table 1 : Electromagnetic Compatibility Information - Electromagnetic emissions**

The device is intended to be used in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.	
Emissions Test	Compliance
Not applicable	Not applicable
CISPR11 Radiated emission	Group 1 Class B
Harmonic emissions CEI 61000-3-2	Not applicable
Voltage fluctuations/flicker emissions CEI 61000-3-3	Not applicable

**Table 2 : Electromagnetic Compatibility Information - Electromagnetic emissions**

The device is intended to be used in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Immunity Test	CEI 60601 Test level	Compliance level
Electrostatic discharge (ESD) CEI 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	contact +/- 8 kV air +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV
Electrostatic transient/burst CEI 61000-4-4	+/- 2 kV for power supply lines 100 kHz repetition frequency +/- 1 kV for input/output lines	N/A
Surge CEI 61000-4-5	+/- 0,5 kV, +/- 1 kV differential mode line-line	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines CEI 61000-4-11	0% UT (100% dip in UT) for 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30% dip in UT) for 25/30 cycles at 0° 0% UT (100% dip in UT) for 250/300 cycles at 0°	N/A
Power frequency (50/60 Hz) magnetic field CEI 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz
<b>NOTE</b> : UT is the A.C. mains voltage prior to application of the test level		

**Table 3 : Electromagnetic Compatibility Information - Electromagnetic emissions**

The device is intended to be used in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Immunity test	CEI 60601 Test level	Compliance level
Conducted RF CEI 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands	N/A
Radiated RF CEI 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m4
Proximity magnetic fields CEI 61000-4-39	30 kHz, 8 A/m, CW 134,2 kHz, 65 A/m Pulse modulation 2,1 kHz 13,56 kHz, 7,5 A/m, Pulse modulation 50 kHz	30 kHz, 8 A/m, CW 134,2 kHz, 65 A/m Pulse modulation 2,1 kHz ; 13,56 kHz, 7,5 A/m, Pulse modulation 50 kHz

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and by reflections from structures, objects and people.

- a. The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are the following: 6 765 MHz to 6 795 MHz ; 13 553 MHz to 13 567 MHz ; 26 957 MHz to 27 283 MHz ; and 40,66 MHz to 40,70 MHz.  
The amateur radio bands between 0,15 MHz et 80 MHz are the following: 1,8 MHz to 2,0 MHz ; 3,5 MHz to 4,0 MHz ; 5,3 MHz to 5,4 MHz ; 7 MHz to 7,3 MHz ; 10,1 MHz to 10,15 MHz ; 14 MHz to 14,2 MHz ; 18,07 MHz to 18,17 MHz ; 21,0 MHz to 21,4 MHz ; 24,89 MHz to 24,99 MHz ; 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- b. The compliance levels in ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that portable/mobile communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in this frequency ranges.
- c. Field strenghts from fixed transmitters, such as base station for radio (cellular/cordless) telephone and land mobile radios, amateur radios, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed transmitters, an electromagnetic site survey should be considered. If the measured field strenght in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary such as re-orienting or relocating the device.
- d. Over the frequency range 150 kHz to 80 MHz, field strenght should be less than 3 V/m.

**Table 4 : Recommended separation distances between RF wireless communication equipments**

The device is intended to be used in the electromagnetic environment in which radiated RF are controlled. The customer or the user of the device can help prevent electromagnetic interferences by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.				
Frequency MHz	Maximum Power W	Distance	Test level CEI 60601	Compliance level
385	1,8	0,3	27	27
450	2	0,3	28	28
710	0,2	0,3	9	9
745				
780				
810	2	0,3	28	28
870				
930				
1 720	2	0,3	28	28
1 845				
1 970				
2 450	2	0,3	28	28
5 240	0,2	0,3	9	9
5 500				
5 782				
<b>NOTE 1</b> : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and by reflections from structures, objects and people.				

### WARNINGS !

- ➔ This device should not be used in the vicinity or on the top of other electronic equipment, such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- ➔ The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement part for internal components, may result in increased emissions or decreased immunity of the equipment or system.

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**SPENGLER SAS**
**Shanghai International Holding Corp. GmbH (Europe)**