LUXAMED®

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Email: info@luxamed.de * www.luxamed.de

- ENGLISH -

Before using the LUXAMED® LED examination light and accessories, read these instructions for use carefully and keep them for later reading.

LUXAMED® LED examination light

Packing:

1x LED examination light

1 x Base

1 x Tripod tube

5 x Foot tubes

1x metal plate, 1x screw

Optional accessories (not included)

Wall holding bracket (Art. Nr. U5.232.014)

Table clamp (Art. Nr. U5.242.014)

Standard rail claw (Art. Nr. U5.252.014)

Purpose:

The LUXAMED® LED examination light is intended exclusively for operation in medical practice. The examination light is operated with the local power grid. It serves for the general and supportive illumination of the examination site. The examination light may only be used by medically trained personnel.

Contraindications:

Dangerous optical radiation may be emitted by this luminaire. Do not look into the lamp during operation. Can be harmful to the eyes.

Warnings:

Please note that the head of the examination light heats up after a longer period of use.

The examination light may only be used if a possible failure does not cause a hazard. Another main lighting must be present.

The examination light must not be used in the MRI area or areas with increased magnetic field and in potentially explosive atmospheres.

Never cover the head of the examination light.

Near flammable anesthetics (class AP) or flammable anesthetics with oxidizing agents (class APG) the luminaire must not be used.

Make sure that the power cord does not bend or become too tense. In addition, care must be taken to ensure that the power cord does not cause a tripping hazard.

Make sure that the examination light is not contaminated during use.

When switching on the luminaire, a one-time short flash may occur due to the integrated current monitoring.

To avoid danger to patients, do not carry out maintenance work when the light is in operation.

The power plug is the disconnecter and must always be easily accessible.

The product is only approved for indoor use, belongs to surge class II and may be used at a maximum altitude of 3000 m.


Do not bend the metal hose too much, never fall below a bending radius of 60 mm (80 mm with model focus).

Any serious incidents involving the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Mounting of the tripod:

First, mount the tripod by inserting the tripod tube into the plastic footboard (1). Then place the metal plate with the elevations in the direction of the tripod tube and attach the screw (2)... Tighten the screw only so far that the foot tubes can still be plugged in laterally (3). The elevations must protrude into the holes in the foot pipes. After attaching the 5 foot tubes, tighten the screw firmly until the tripod is stable.

Now press the arm of the examination light with the recording into the tripod. The luminaire must be fixed to the tripod base by the lateral pin screw. To do this, tighten the pin screw clockwise with the supplied Allen wrench.

The examination light can now be plugged into the local power grid. By pressing the main switch, the luminaire can be switched on or off.

Make sure that a minimum distance to the inspection site of 200 mm is maintained.

Focusing on the model Focus: By turning the black plastic screw, the light spot can be varied in size and illuminance.

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Desinfektion:

Turn off the light and then unplug it before each cleaning.

Cleaning with a soft and lint-free cloth, moistened with a little alcohol if necessary. Make sure that you do not touch any electrical contacts. The luminaire must be cleaned after each application.

Disinfection:

According to the recommendation of the Commission for Hospital Hygiene and Infection Prevention (KRINKO) of the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM):

Only disinfectants with proven antibacterial, antiviral and fungicidal efficacy should be used. Bundesgesundheitsblatt 2012 · 55:1244-1310; DOI 10.1007/s00103-012-1548-6
In principle, CE marked alcohol-based instrument disinfectants (2-propanol) are suitable.

Please follow the instructions for use of the disinfectant used.

Wipe the lamp head with a disinfectant-soaked cloth. Observe the exposure time according to the instructions for use of the disinfectant. Then wipe away any residue

Sterilization:

Sterilization is not permitted.

Maintenance:

No maintenance is required. In case of defect, the luminaire can be repaired by the manufacturer.

Control and functional testing:

Check the integrity of the luminaire and power cord before each application. By switching it on, you can see if the examination light is working correctly. When switched on, there may be a short flare-up due to current peaks, which are filtered by the power supply, which does not affect the function and safety. This protects the LEDs and increases the service life of the LED bulb.

Disposal:

The examination light with arm must be fed to a separate collection of electrical and electronic equipment. The tripod can be fed to the scrap metal.

Technical data:

Weight (incl. tripod):
Mains connection:
Lamp life:
Product lifetime:
Illuminance:
Color temperature:
Min. CRI
Maximum value total irradiation

6.2 kg
100–240 V AC, 50/60 Hz, 0.6 A
approx. 50,000 hours
min. 5 years
~ 50,000 lux at 300 mm distance
approx. 4,000 K

LED examination light focus

6.6 kg
100–240 V AC, 50/60 Hz, 0.6 A
approx. 50,000 hours
min. 5 years
~ 80,000 lux at 300 mm distance
approx. 5,000 K
92.4
Ee = 1.4e+00 W/m²

Environmental conditions

Temperature:
Rel. Humidity:
Air pressure:

+ 5°C to + 40°C
10% to 90%
700 hPa to 1070 hPa

Storage/transport conditions:

Temperature
Rel. Humidity:
Working height max.:
IP protection class:
Length of flexible arm:

-15° to +45°C
10% to 90%
1.650mm
IP 20
900mm


Minimum bending
radius:

60mm

80 mm

Symbol explanation

	Manufacturer
	Date of manufacture
	Follow the instructions for use
	Warnings
	Separate collection of electrical and electronic equipment (European WEEE Directive)
	Compliant with EU regulations
	The device is a medical device
	Interference is possible in the vicinity of devices bearing the figurative sign.
	Protection class II
	may only be used in dry rooms
	Compliant with UK regulations
	Swiss authorized representative for medical devices

Warranty conditions

For this product we provide 2 years factory warranty from the date of purchase. The date of purchase must be proven by invoice. Within the warranty period, defects resulting from material or manufacturing defects will be remedied free of charge. The warranty service does not extend the warranty period to the entire device, but only to the replaced components. Excluded from the warranty is wear and tear due to wear and tear, transport damage and all damage caused by improper handling (e.B. non-observance of the instructions for use) or due to interventions by unauthorized persons. The guarantee does not establish any claims for damages against us. The buyer's statutory claims for defects in accordance with § 437 BGB (German Civil Code) are not restricted. In case of claiming the warranty, the device must be sent to the manufacturer together with the original proof of purchase :

LUXAMED GmbH & Co. KG
Daniel-Weil-Straße 3
D-89143 Blaubeuren
GERMANY

After-sales service: Warranty and repair work must be carried out by trained and authorized personnel. The device must not be modified without the manufacturer's permission. Send the device carefully packed and sufficiently franked to your authorized dealer or directly to:

LUXAMED GmbH & Co. KG
Abteilung Service
Daniel-Weil-Straße 3
D-89143 Blaubeuren
GERMANY
Phone: + 49 (0) 7344 92905-0

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EMC Guidelines

Guidance and manufacturer' declaration - electromagnetic emissions - for LUXAMED® examination lamp.

Emission

Test	Limit	Electromagnetic environment - guidance
Conducted emission	CISPR 11, Group 1, Class B	Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated emission	CISPR 11, Group 1, Class B	
Harmonic current emissions	IEC 61000-3-2, Class A	Device is suitable for direct connection to the public low-voltage power supply network.

Immunity test levels

Test	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: ±8 kV Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM field (IEC 61000-4-3)	80-2700 MHz; 1kHz AM 80 %; 10 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ for 800 MHz to 2,7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>

<p>Proximity fields from RF wireless communications equipment (IEC 61000-4-3)</p>	<p>385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM + 5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance 30 cm.</p>
---	---	--

<p>Electrical fast transients / bursts (IEC 61000-4-4)</p>	<p>Power lines: 2 kV; 100 kHz repetition frequency Signal lines: 1 kV; 100 kHz repetition frequency</p>	<p>Mains power quality should be that of a typical environment.</p>
<p>Surges (IEC 61000-4-5)</p>	<p>L-N: 1kV at 0°, 90°, 180°, 270° L-PE, N-PE: 2 kV at 0°, 90°, 180°, 270°</p>	<p>Mains power quality should be that of a typical environment.</p>
<p>Conducted disturbances induced by RF fields (IEC 61000-4-6)</p>	<p>0.15-80 MHz; 1kHz AM 80 %; 3 Vrms, 6 Vrms in ISM and amateur radio band</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2\sqrt{P}$ for 150 kHz to 80MHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>
<p>Rated power frequency magnetic fields (IEC 61000-4-8)</p>	<p>30 A/m, 50 Hz</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

<p>Voltage dips / Voltage interruptions (IEC 61000-4-11)</p>	<p>0 % UT for 0.5 cycle at 0° ,45°,90°,135°,180°,225°,270°,315° 0 % UT for 1 cycle at 0° 70 % UT for 25/30 cycles at 0° 0 % UT for 250/300 cycles at 0°</p>	<p>Mains power quality should be that of a typical environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or battery.</p>
--	--	--

www.luxamed.de

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LUXAMED GmbH & Co. KG. KG.