

CB48

HEATING BLANKETS FOR OPERATING AND RECOVERY ROOM



**MADE IN
FRANCE**



- Reusable!
- Easy to use
- Safe
- Made In France

Clinical studies demonstrated that perioperative hypothermia may cause complications for patients. CB48 heating blankets are patient active warming systems dedicated to operative rooms, post-interventional monitoring rooms (SSPI) and for warming up patients prior to surgery.

CE
0459

GSH

30 rue Jean de Guiramand
13290 Aix en Provence | FRANCE
www.gsh-med.fr

Contact :
+33 (0) 4 42 90 31 31 - info@gsh-med.fr

CB48 blankets are not using pulsed air technology but an integrated circuit in a cartridge (prevention of noises and particles). Costs are controlled as they are reusable and do not need consumables. In this way, they prevent hypothermia in patients undergoing surgery.

Blankets are controlled by a BCB48 electronic control unit with LCD display.



SUSTAINABLE

- Adapted technology: integrated circuit in a cartridge, no pulsed air, preventing noise and presence of particles in the operating environment
- Inertia when power is cut off, keeping the patient warm for some time
- Reusable: reduces consumable waste and is more economical



EASY TO USE:

- Reusable
- Various sizes to suit every patient.

	<i>Width</i>	<i>Length</i>
CB48-2 +BCB48: Half body/pediatric	80cm	100cm
CB48-3 +BCB48: Upper body	50cm	150cm
CB48-11 +BCB48: Full body	100cm	150cm

- Select temperature from 30 to 40°C with the BCB48 knob and switch on the machine.
- The blankets are made with PU cover, which is very easy to clean and disinfect after use.



SAFE:

- Precision and safety with 8 independent temperature sensors.
- Equipped with visual and audible alarms to call nurses in case of:
 - Heating element or blanket conductor breaking
 - Temperature sensor default
 - Temperature variation too high
 - Temperature above 43°C
 - Power outage
 - Main board default
- Low voltage direct current power: 48V, to prevent from monitors interferences.



MADE IN FRANCE:

We develop and produce patient warming systems since 1991 in our production site of Saint Leonard de Noblat, near Limoges - FRANCE !

References

DESCRIPTION	REFERENCE
Heating Blanket - Half body/Pediatric	CB48-2+BCB48
Heating Blanket - Adult (half body)	CB48-3+BCB48
Heating Blanket - Adult (full body)	CB48-11+BCB48



Accessories:

DESCRIPTION	REFERENCE
Cover holder for serum stand	SUPPS



ECM FRENCH QUALITY

CE MARKING 0459

Warming cots are class IIb medical devices. They comply with current European regulations.

CUSTOMER TECHNICAL SUPPORT CENTRE

GSH provides maintenance and repair services for heating blankets. Contact our after-sales service on : +33 5 55 56 32 40.

ECM WARRANTY

ECM heating blankets are guaranteed for 12 months (under normal conditions of use). The warranty does not cover maintenance faults made by the user.