



Spengler

SINCE 1907



# *Cardio prestige II*

SINGLE CHESTPIECE  
PROFESSIONAL STETHOSCOPE

CE



Spengler

Spengler S.A.S | ZI La Limoise | Rue Robert Maréchal | 36100 ISSOUDUN | FRANCE

[www.spengler.fr](http://www.spengler.fr)

Tel. : +33 (0)4 42 52 89 80

Mail : [export@spengler.fr](mailto:export@spengler.fr)

# PROFESSIONAL STETHOSCOPE

## Cardio prestige II SINGLE CHESTPIECE

### *The excellence.*

*With its heavy chestpiece equipped with 'premium acoustic' technology, the new Cardio prestige II stethoscope provides the best results.*



#### PREMIUM ACOUSTICS

Precision-machined, the new design for the CARDIO PRESTIGE II receiver is inspired by musical instruments such as the French horn. Its flared shape allows sound waves to be captured better.



#### STAINLESS STEEL CHESTPIECE

Stainless steel offers the best acoustic properties, so it was selected for use in the chestpieces for the CARDIO PRESTIGE II stethoscopes. In addition, we increased the thickness where necessary in order to refine the acoustics.



#### ADAPTED TO CUFFS

Adapted to the cuffs, the CARDIO PRESTIGE II stethoscopes are designed with thin edges, perfectly matching the ends of the tensiometer cuff.



#### FLOATING MEMBRANE

Present on both chestpiece:

- Light pressure to hear low frequencies
- Firm pressure to hear high frequencies



#### OVERSIZE TUBE

The CARDIO PRESTIGE II is equipped with a wide double-duct tubing and is perfectly insulated against noise.



#### DESIGNATION

● Cardio Prestige® II stethoscope single chestpiece

#### REFERENCE

506 010

### SPENGLER QUALITY

#### CE mark

European directive 93/42/CEE : stethoscopes: class I. The Spengler stethoscopes comply with the stipulations of the guidelines for safety, health, hygiene and protection of the environment.

#### After-sales service

Spengler carries out the maintenance and repairs.

#### Legal guarantee

According to French law, the Spengler stethoscopes are guaranteed for life against any hidden manufacturing defect.

#### Spengler guarantee

Spengler guarantees the stethoscope Cardio Prestige II for 7 years from date of purchase. The guarantee does not apply to damage caused by improper

handling, accidents or alterations made to the instrument by third parties. Spengler guarantees the ear tips, diaphragm and non-chill ring for 1 year from date of purchase.

The guarantee does not apply to fair wear and tear of the product/damages and breakdowns due to shock, misuse or neglect/any product disassembled, repaired or modified by the user or unauthorised persons.



[www.spengler.fr](http://www.spengler.fr)

Spengler S.A.S | ZI La Limoise | Rue Robert Maréchal | 36100 ISSOUDUN | FRANCE

Tel. : +33 (0)4 42 52 89 80  
Mail : [export@spengler.fr](mailto:export@spengler.fr)