

## PNEUMATIC TOURNIQUET CUFFS

Arm cuffs and lower limb cuffs are used with pneumatic tourniquets for orthopaedic surgery.

**DDM tourniquet cuffs are made in France!**



CE

1/3

### AN ADAPTED CUFF RANGE FOR ALL TYPES OF PATIENTS

#### Monobloc PU arm and lower limb cuffs

Monobloc PU arm and lower limb cuffs, decontaminable with 1 or 2 bladders, made of a single Monobloc unit with a ribbon to secure the Velcro or a fastening strap (optional) and featuring a range of sizes identifiable by color-coding of the tubing.

Available in straight and conical design to better fit the limb.



	WITH STRAP		WITH STRAP		WITH STRAP	
	1 BLADDER - STRAIGHT		1 BLADDER - CONICAL		2 BLADDERS - STRAIGHT	
Infant arm cuff	GBDM101	GBDSM101	/	/	/	/
Child arm cuff	GBDM102	GBDMS102	/	/	GBDM202	GBDMS202
Small adult arm cuff	GBDM103	GBDMS103	/	/	/	/
Adult arm cuff	GBDM104	GBDMS104			GBDM204	GBDMS204
Large adult arm cuff	/	/	GBCM105	GBCMS105	/	/
					2 BLADDERS - CONICAL	
Small adult lower limb cuff	GCDM103	GCDMS103	/	/	/	/
Adult lower limb cuff	GCDM104	GCDMS104	GCCM104	GCCMS104	GCCM204	GCCMS204
Large adult lower limb cuff	GCDM105	GCDMS105	GCCM105	GCCMS105	/	/
XL adult lower limb cuff	GCDM106	GCDMS106	GCCM106	GCCMS106	/	/
XXL adult lower limb cuff	GCDM107	GCDMS107	GCCM107	GCCMS107	/	/

#### Autoclavable polyester arm and lower limb cuffs

Disinfectable and autoclavable polyester arm and lower limb cuffs are made of a removable silicone pouch for easy autoclaving.










Available in straight and conical design to better fit the limb.



	1 BLADDER - STRAIGHT	2 BLADDERS - STRAIGHT	LOWER LIMB CUFFS - CONICAL
Infant arm/lower limb cuff	GBS101	/	/
Child arm cuff	GBS102	GBS202	/
Small adult arm cuff	GBS103	/	/
Adult arm cuff	GBS104	GBS204	/
Child lower limb cuff	GCS102	/	/
Small adult lower limb cuff	GCS103	/	/
Adult lower limb cuff	GCS104	GCS204	GCCS104
Large adult lower limb cuff	GCS105	/	GCCS105
XL adult lower limb cuff	GCS106	/	GCCS106

### Connectors

#### CLIP CONNECTION SYSTEM WITH O-RING

REFERENCE	IMAGE	MATERIAL	DESIGNATION
A11012M		METAL	4.8 mm fluted in-line female coupler
A11016		ACETAL	4.8 mm fluted in-line female end-piece with sealing valve
A11008M		METAL EPDM	4.8mm metal male connector with O-ring
A11010		ACETAL EPDM	4.8 mm fluted in-line male end-piece with sealing valve
A11322		EPDM	O-ring for male connector
A20001		METAL	3 m blue extension cable + male and female connectors
A20006		METAL	3 m red extension cable + male and female connectors
A20003		SILICONE + ACETAL	3 m blue autoclavable, underwater extension cable + male and female connectors
A20005		SILICONE + ACETAL	3 m red autoclavable, underwater extension cable + male and female connectors

### CLEANING & DISINFECTION

Monobloc arm and lower limb cuffs cannot be autoclaved. They can be dry cleaned or cleaned upon immersion taking care to maintain the connector (on the tube) out of the water.

Autoclavable cuffs can go through autoclaving according to standard protocols.

#### DDM FRENCH QUALITY

##### CE mark

Tourniquet cuffs are Class I medical devices. They comply with current European regulations.

##### After-Sales Service

GSH provides tourniquet cuffs maintenance and repair services. Contact our After-Sales Service on: +33 4 42 90 31 31.

##### Warranty

Monobloc cuffs are guaranteed for 2 years, and autoclavable cuffs are guaranteed for 1 year (under normal conditions of use). The warranty does not cover faulty maintenance by the user.