

# CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM

### Ref. DUCRYO03







**IEDICAL DEVICI** 

CE



29/08/2024 19:48:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand 13290 Aix-en-Provence | France www.gsh-med.fr **↓** +33(0)4 42 90 31 31**⋈** info@gsh-med.fr

(1/2



## CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM

# Utermöhlen Since 1880

2/2

#### Ref. DUCRYO03

Effectively treats benign skin lesions

### **Description**

- Cryotherapeutic agent, allowing to deliver a temperature of -

55  $^\circ$  C in less than 15 seconds

Application time between 15 and 40 seconds

Complete set including:

A 170 ml aerosol can

2mm and / or 5 mm applicator

sticks, usable depending on the different lesions

An instruction manual

- The sale of this medical equipment is reserved for health institutions, doctors and doctors' offices

## **Features**

CE Marking	CE0344	Medical device	Class lla
Barcode	8717484005835	HS Code	90189084
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	14
Packed unit / Length (cm)	14	Packed unit / Weight	0.3
SPECIFIC STORAGE CONDITIONS			
Maxtemperature	25		

OTHER FEATURES

29/08/2024 19:48:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand 13290 Aix-en-Provence | France www.gsh-med.fr +33(0)4 42 90 31 31info@gsh-med.fr