

## FOOT STEP - 1 STEP - EPOXY - HOLTEX



Ref. OS0051E



29/08/2024 18:32:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## FOOT STEP - 1 STEP - EPOXY - HOLTEX



Ref. OS0051E

### Description

- Structure in epoxy lacquered steel tube
- Walk covered with non-slip plastic
- Maximum load: 150 kg
- Dimensions: (L)40 x (W)24.5 x (H)22 cm

Stainless steel or Epoxy: what difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168098	HS Code	9402900000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	31	Bare Unit / Width (cm)	24.5
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	22
Packed unit / Length (cm)	42	Bare Unit / Length (cm)	40
Packed unit / Weight	4	Bare Unit / Weight (kg)	3.4
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 29/08/2024 18:32:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)