


## LITHIUM REPLACEMENT BATTERY FOR DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. SSACCE010



30/08/2024 15:45:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LITHIUM REPLACEMENT BATTERY FOR DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. **SSACCE010**

### Description

The battery can be stored for 8 years in its original packaging while maintaining its original performance. It is not necessary to replace it after use on patient if the device does not indicate it.

- Compatibility: FRED PA-1 defibrillator

- Technology: Li/MnO<sub>2</sub>- Power: 15

V/2.8 Ah

- Battery life: more than 140 shocks at a maximum energy level and 6 hours of operation, if the optimal storage / use temperatures of the device are between 15 - 25 °C- Standby life: • 6 years with weekly self-tests for the FRED PA-1

• 3 years with weekly self-tests for the FRED®

© PA-1 Online

- Storage and transport temperature: 5 to 35°C (48h max between -20° to 5°C and 35° to 60°C)

### Features

CE Marking	NONE	Medical device	No
Bar code	3661041009222	HS Code	90189010
Warranty (months)	72	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	6	Packed unit / Height (cm)	3
Packed unit / Length (cm)	19	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.28		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 15:45:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)