

## STRIPS FOR SPENGLER TB100 GLUCOMETER (50UNITS)



Ref. 745010



CE  
0537

DISPOSITIF MEDICAL  
**IVD**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 YEARS  
GUARANTEE

STERILE

 30/08/2024 20:29:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## STRIPS FOR SPENGLER TB100 GLUCOMETER (50UNITS)



Ref. 745010

### Description

- Bottle of 50 strips for Glucometer SPENGLER TB 100
- Strips to be used in addition to the Complete Kit ref. 745000.
- For in vitro diagnosis only
- Blood glucose test on capillary blood
- Use within 3 months after opening (store between 4 ° C and 30 ° C)The TB 100

test strips are used with the Spengler TB 100 glucometer for the quantitative analysis of the level of glycemia (sugar) in whole blood which is taken at the fingertip, in the palm or forearm. TB 100 test strips are designed for out-of-body analysis (in vitro diagnostics). The blood glucose monitoring system, the TB 100 blood glucose meter, is intended for use at home and in the workplace to monitor blood glucose levels. Remember to throw all your care items (strips, cotton and cotton wool) in the collection boxes for infectious waste.

Refund BVG code: 6111350

### Features

CE Marking	CE0537	Medical device	DMDIV (diagnostic In Vitro)
Bar code	3700446041656	HS Code	90183190
Warranty (months)	24	Sterile	Yes
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	4	Packed unit / Height (cm)	4
Packed unit / Length (cm)	7	Packed unit / Weight	0.035
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	4	Max temperature	30
<b>OTHER FEATURES</b>			
Code LPP	6111350	Packed unit / Quantity	1

30/08/2024 20:29:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)