

GREY FUSE[®], CLEAR FUSE[®] & EASY FUSE[®] PRESSURE INFUSION CUFFS

All our pressure infusion cuffs are manufactured in France, in our Manufacture in Castelculier (near Agen)!



GREY FUSE[®]



CLEAR FUSE[®]



EASY FUSE[®]

- French quality.
- Optimized monitoring
- Reliability



GREY FUSE[®], CLEAR FUSE[®] & EASY FUSE[®]

GREY FUSE[®], CLEAR FUSE[®] AND EASY FUSE[®] pressure infusion cuffs are medical devices designed for daily use in operating theaters, intensive care units or emergency transport vehicles. They can be used to **manually accelerate infusion** of blood and solutions when gravity alone cannot achieve the desired flow rates.

GREY FUSE[®], CLEAR FUSE[®] AND EASY FUSE[®] pressure infusion cuffs come in 3 sizes: 500cc, 1000cc, 3000cc, suitable for adapting to different sizes of solution or blood bags.



French quality:

All our pressure infusion cuffs are manufactured in France in our Manufacture in Castelculier (near Agen)!




Optimized monitoring:


- Visualization of the solution bag:
 - thanks to a transparent fabric composed of polyamide weave and mesh for GREY FUSE[®] and EASY FUSE[®]
 - thanks to transparent polyurethane material that wraps around the CLEAR FUSE[®] bag
- Pressure gauge accuracy of ± 3 mmHg for precise pressure monitoring for the GREY FUSE[®] and CLEAR FUSE[®], pressure indicator with automatic overpressure regulator on EASY FUSE[®]
- Red zone on the dial at 300 mmHg to prevent overpressure.
- Controlled decompression facilitated by a single-way valve for GREY FUSE[®] and CLEAR FUSE[®], and a 3-way valve for EASY FUSE[®].




Reliability:

- These robustly-designed cuffs provide great resistance as they are made of TPU-coated fabric (GREY FUSE[®] and EASY FUSE[®]) or high-frequency-welded TPU membrane (CLEAR FUSE[®]).
- Suitable for all uses:
 - compatible with all blood/solution bags (500cc, 1000cc, 3000cc)
 - Reusable (GREY FUSE[®] and CLEAR FUSE[®]) and single-use (EASY FUSE[®])
- Very easy to clean or disinfect thanks to easy disassembly of the different parts (only for reusable cuffs).

<p>GREY FUSE[®]</p> 	LATEX-FREE REUSABLE NON-STERILE 3 SIZES		
	REFERENCE	SIZE	BAG DIMENSIONS IN MM
	M20500	500 cc	250 x 162
	M21000	1000 cc	365 x 170
	M25000	3000 cc	561 x 350

<p>CLEAR FUSE[®]</p> 	LATEX-FREE REUSABLE NON-STERILE 3 SIZES		
	REFERENCE	SIZE	BAG DIMENSIONS IN MM
	M30500	500 cc	220 x 191
	M31000	1000 cc	333 x 222
	M34000	3000 cc	420 x 440

<p>EASY FUSE[®]</p> 	LATEX-FREE REUSABLE NON-STERILE 3 SIZES		
	REFERENCE	SIZE	BAG DIMENSIONS IN MM
	M10500	500 cc	250 x 162
	M11000	1000 cc	365 x 170
	M13000	3000 cc	338 x 235

FRENCH QUALITY DDM

CE mark

Pressure infusion cuffs are Class I medical devices. They comply with current European regulations.

Customer Technical Support Centre

GSH provides pressure infusion cuffs maintenance and repair services. Contact our After-Sales Service on: +33 4 42 90 31 31.

Warranty

DDM pressure infusion cuffs are guaranteed for 24 months (under normal conditions of use). The warranty does not cover faulty maintenance by the user.