

CLEAR FUSE® PRESSURE INFUSION CUFFS

TECHNICAL SHEET


All our pressure infusion cuffs are manufactured in France, in our Manufacture in Castelculier (near Agen)!



The special design of this pressure cuff allows it to be wrapped around the liquid bag, making it very easy to use and clean. Pressure is controlled by a pressure gauge. Vacuum is controlled by a tap.

CE

1/2



Letter	Reference	Description	Feature
A	MS3050	Cuff alone 500 cc	Polyurethane
A	MS3100	Cuff alone 1000 cc	Polyurethane
A	MS34000	Cuff alone 3000 cc	Polyurethane
B	A10807	Pressure gauge	ABS + plexiglas
C	PM50B	Pressure gauge protector	Silicon
D	A10259	Bulb tap	Aluminium + brass
E	A10219	Bulb	PVC
D/E	A10226	Complete Bulb	ABS + PVC

CLEANING & DISINFECTION	Surface cleaning	Whole device
	Disinfection with decontamination solution	Whole device
STORAGE CONDITIONS	Packaging type	Origin
	Local	Ventilated
	Temperature	-10° to 40 °C
	Hygrometry	30 to 40 %
	Atmospheric pressure	500 to 1060 hPa

Description	Size	Bladder dimensions (mm)		Reference
		Lenght	Width	
<i>Transparent pressure infusion cuff</i>	500 cc	220	191	M30500
	1 000 cc	333	222	M31000
	3 000 cc	420	440	M34000

FRENCH QUALITY DDM

CE mark

Pressure infusion cuffs are Class I medical devices. They comply with current European regulations.

Customer Technical Support Centre

GSH provides pressure infusion cuffs maintenance and repair services. Contact our After-Sales Service on: +33 4 42 90 31 31.

Warranty

DDM pressure infusion cuffs are guaranteed for 24 months (under normal conditions of use). The warranty does not cover faulty maintenance by the user.