

GREY FUSE® PRESSURE INFUSION CUFFS

TECHNICAL SHEET

All our pressure infusion cuffs are manufactured in France, in our Manufacture in Castelculier (near Agen)!






Grey Fuse pressure cuffs are made from a single unit:

- robust design
- high-quality technical fabric
- easy cleaning and decontamination
- wrist strap for attaching the device to the serum stand on the M25000.
- hook to hold the solute bag on the M25000 due to the weight of the solute bags
- transparent polyamide weave and mesh fabric for easy viewing of the solute bag

((

1/2

PRESSURE INFUSION CUFFS



Letter	Reference	Description	Feature		
G	M20511	Cuff alone 500 cc	PU-coated fabric + transparent polyamide fabric		
G	M20011	Cuff alone 1000 cc	PU-coated fabric + transparent polyamide fabric		
G	M25011	Cuff alone 3000 cc	PU-coated fabric + transparent polyamide fabric		
E	A10808	Pressure gauge	ABS + PC glass + TPU protector		
С	A10259M	Bulb tap	Aluminum + chrome-plated brass		
A-B-C	A10226M	Complete bulb with metal tap	Phthalate-free ABS + PVC		
D	-	Tubing	Phthalate-free PVC		
F	-	Transparent mesh	Polyamid		
Н	-	Strap	Polyester		
I	-	Hand Strap	Polyester		

CLEANING &	Surface cleaning	Whole device
DISINFECTION	Disinfection with decontamination solution	Whole device
	Packaging type	Origin
CTOD 4 CT	Local	Ventilated
STORAGE CONDITIONS	Temperature	-10° to 40 °C
CONDITIONS	Hygrometry	30 to 40 %
	Atmospheric pressure	500 to 1060 hPa

Description	Size		dimensions nm)	Reference
		Lenght	Width	
	500 cc	250	162	M20500
Monobloc pressure infusion cuff	1 000 cc	365	170	M21000
iniusion cun	3000 cc	561	350	M25000

FRENCH QUALITY DDM

CE mark

Pressure infusion cuffs are Class I medical devices. They comply with current European regulations.

Customer Technical Support Centre

GSH provides pressure infusion cuffs maintenance and repair services. Contact our After-Sales Service on: +33 4 42 90 31 31.

Warranty

DDM pressure infusion cuffs are guaranteed for 24 months (under normal conditions of use). The warranty does not cover faulty maintenance by the user.

FT_GREYFUSE_IVB_20240411_EN