

## WALL STAND FOR DEFIBRILLATOR FRED PA-1 SCHILLER




Ref. SSACCE002



 30/08/2024 15:40:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## WALL STAND FOR DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. **SSACCE002**

### Description

The docking station is thermoformed to fit the curves of the FRED PA-1. It allows a grip of the device. Discreet, it blends into all interiors.

Compatibility: FRED PA-1, FRED PA-1 Online and FRED PA-1 Trainer defibrillators - Indoor installation only.  
- 2 slides to fix it to the wall.

### Features

CE Marking	NONE	Medical device	No
Bar code	3661041009147	HS Code	90189010
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	16	Bare Unit / Width (cm)	16.5
Packed unit / Height (cm)	13	Bare Unit / Height (cm)	11.5
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	27
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.25
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 15:40:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)