



SMARTLED
SPENGLER X LUXAMED EDITION®
The innovative pocket size otoscope MADE IN GERMANY



CE



Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr

Tel. : +33 (0)4 42 52 89 89
Mail : contact@spengler.fr

SMARTLED SPENGLER X LUXAMED EDITION®



LIGHT ILLUMINANCE
10 000 Lux



INNOVATION LED
6 Rings of LEDS (concentric) placed at the front of the otoscope allowing an unrestricted field of view.



LONG LIFESPAN
About 100,000 hours with constant current



Brightness dimmer
A long press is all it takes to adjust the brightness.



AUTOMATIC STOP
After 3 min of no use



ZOOM X3
With a swivel magnifier with magnification x3.

TECHNICAL CHARACTERISTICS

Head	Insufflation port	For pneumatic otoscopy
	Pivoting magnifying glass	3X magnification Material : PCNT
Handle	Material	Aluminium and reinforced fiberglass
	Technology	LED
Éclairage	Voltage	2,5 Volts
	Lifetime	100 000 hours
Colour temperatures	4.200 k	
CRI	>90	
Illuminating	approximately 10 000 LUX	
Power	Batteries	2 alkaline AAA (supplied)
Accessories	Supplied	20 grey disposable specula including 10 of 2.5 mm and 10 of 4.0mm
Total length	16,2 cm	
Storage box	Dimensions	24.5 cm x 11.3 cm x 4.1 cm
SPENGLER GUARANTEE	3 YEARS	Head and handle (under normal use conditions), Excluding consumables (specula, batteries). Against any hidden manufacture defect.

DESCRIPTION	CODE
Otoscope SmartLed Spengler x LUXAMED edition - Black	A1.416.112
Otoscope SmartLed Spengler x LUXAMED edition - Bronze	A1.416.184
Otoscope SmartLed Spengler x LUXAMED edition - Blueberry	A1.416.194
Otoscope SmartLed Spengler x LUXAMED edition - Pediatric	A1.516.514

Accessories

Auricular speculum 4 mm (box of 250)	L20848
Auricular speculum 2.5 mm (box of 250)	L20747



SPENGLER QUALITY

Marquage CE

The Spengler otoscopes including the SmartLED Spengler x Luxamed edition are medical devices of class I. They comply with the European Regulation 2017/745, perfectly meeting the general safety and performance requirements.

Customer service department

Spengler SAS provides maintenance and repair of equipment. Head over to the official website in the Contact section.

Spengler 3 year warranty

(under normal conditions of use): Free repair or replacement of parts found defective (defect of operation due to a material, manufacturing or design defect) by our after-sales service. The warranty does not cover maintenance defects on the part of the user.



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guirmand | 13290 Aix en Provence | FRANCE

Tel. : **+33 (0)4 42 52 89 89**

Mail : contact@spengler.fr