

GUNNAR DRAWING FORCEPS 18.5 CM

Réf. IP54818



GUNNAR DRAWING FORCEPS 18.5 CM



Réf. IP54818

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040028330
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,7	Width packaging (cm):	
Length (cm) :	18,4	Length packaging (cm) :	
Height (cm) :	2,8	Height packaging (cm) :	
Weight (kg) :	0,098	Weight packaging (kg) :	0,098
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:27 - This article is available on our website : www.gsh-med.fr with the following reference IP54818

All characteristics are given for information only and may be changed