

KERRISSON PUNCH FORCEPS UP 90° 20 CM X 4 MM

Réf. IP46841



KERRISSON PUNCH FORCEPS UP 90° 20 CM X 4 MM



Réf. **IP46841**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026411
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,7	Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :	1,4	Height packaging (cm) :	
Weight (kg) :	0,150	Weight packaging (kg) :	0,182
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:16 - This article is available on our website : www.gsh-med.fr with the following reference IP46841

All characteristics are given for information only and may be changed