

FERRIS-SMITH-KERRISSON FORCEPS DOWNWARD CUTTING 90° - 180 MM - 5 MM

Réf. IP46849



FERRIS-SMITH-KERRISSON FORCEPS DOWNWARD CUTTING 90° - 180 MM - 5 MM



Réf. IP46849

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026497
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,5	Width packaging (cm):	12,5
Length (cm) :	27,6	Length packaging (cm) :	27,6
Height (cm) :	0,8	Height packaging (cm) :	0,8
Weight (kg) :	0,236	Weight packaging (kg) :	0,236
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:16 - This article is available on our website : www.gsh-med.fr with the following reference IP46849

All characteristics are given for information only and may be changed