

TYDING FORCEPS 20 CM

Réf. IP20130



TYDING FORCEPS 20 CM



Réf. IP20130

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe Ir | |
| Marquage | CE | Barcode | 3661040023496 |
| Guarantee (month): | 24 | Customs Code | 9018.9090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 7,2 | Width packaging (cm): | |
| Length (cm) : | 21,4 | Length packaging (cm) : | |
| Height (cm) : | 0,3 | Height packaging (cm) : | |
| Weight (kg) : | 0,063 | Weight packaging (kg) : | 0,063 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 17:06 - This article is available on our website : www.gsh-med.fr with the following reference IP20130

All characteristics are given for information only and may be changed