

## BENGOLEA FORCEPS STRAIGHT 1X2 TEETH 24 CM

Réf. IP42024



# BENGOLEA FORCEPS STRAIGHT 1X2 TEETH 24 CM



Réf. IP42024

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion. The steels come into contact with human tissues and must therefore comply with international standards and regulations in force. Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040025384
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,8	Width packaging (cm):	
Length (cm) :	24,0	Length packaging (cm) :	
Height (cm) :	0,4	Height packaging (cm) :	
Weight (kg) :	0,066	Weight packaging (kg) :	0,066
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP42024

All characteristics are given for information only and may be changed