

## ANTIMAGNETIC AA DISSECTING FORCEPS

Réf. IP27114



## ANTIMAGNETIC AA DISSECTING FORCEPS



Réf. IP27114

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe Ir               |               |
| Marquage                   | CE        | Barcode                 | 3661040024424 |
| Guarantee (month):         | 24        | Customs Code            | 9018.9090     |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 3,6       | Width packaging (cm):   |               |
| Length (cm) :              | 16,9      | Length packaging (cm) : |               |
| Height (cm) :              | 0,5       | Height packaging (cm) : |               |
| Weight (kg) :              | 0,020     | Weight packaging (kg) : | 0,022         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



30/03/2022 17:14 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP27114

All characteristics are given for information only and may be changed