

## DOYEN MALLETT SINGLE HEAD 21 CM 500 GR (STAINLESS STEEL)

Réf. IM02350



# DOYEN MALLET SINGLE HEAD 21 CM 500 GR (STAINLESS STEEL)



Réf. IM02350

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040020402
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,6	Width packaging (cm):	12,6
Length (cm) :	24,8	Length packaging (cm) :	24,8
Height (cm) :	5,0	Height packaging (cm) :	5,0
Weight (kg) :	0,628	Weight packaging (kg) :	0,628
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:49 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IM02350

All characteristics are given for information only and may be changed