

HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°29/30

Réf. IS01429



HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°29/30



Réf. IS01429

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029689
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,7	Width packaging (cm):	
Length (cm) :	21,3	Length packaging (cm) :	
Height (cm) :	2,7	Height packaging (cm) :	
Weight (kg) :	0,146	Weight packaging (kg) :	0,146
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:35 - This article is available on our website : www.gsh-med.fr with the following reference IS01429

All characteristics are given for information only and may be changed