

## SEBILEAU RASPATORY 15 CM 3 MM

Réf. IR09103



# SEBILEAU RASPATORY 15 CM 3 MM



Réf. IR09103

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029375
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,5	Width packaging (cm):	1,5
Length (cm) :	15,6	Length packaging (cm) :	15,6
Height (cm) :	0,8	Height packaging (cm) :	0,8
Weight (kg) :	0,037	Weight packaging (kg) :	0,037
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:34 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IR09103

All characteristics are given for information only and may be changed