

3 INCHES DISPOSABLE PADS FOR TDH39 AND DD45 HEADPHONES (X500)

Réf. DSIPROT01



3 INCHES DISPOSABLE PADS FOR TDH39 AND DD45 HEADPHONES (X500)

Réf. **DSIPROT01**

FEATURES			
Medical device	No		
Marquage	CE	Barcode	
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	30,0
Length (cm) :		Length packaging (cm) :	30,0
Height (cm) :		Height packaging (cm) :	15,0
Weight (kg) :		Weight packaging (kg) :	0,800
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



01/07/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference DSIPROT01

All characteristics are given for information only and may be changed