


## S/S BOWL WITH SPOUT DIAM 160 MM - 1 LITER

Réf. AI15160



## S/S BOWL WITH SPOUT DIAM 160 MM - 1 LITER



Réf. AI15160

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019413
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	16,3	Width packaging (cm):	16,3
Length (cm) :	16,3	Length packaging (cm) :	16,3
Height (cm) :	8,0	Height packaging (cm) :	8,0
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:46 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference AI15160

All characteristics are given for information only and may be changed