


SCALE RANGE 300 KGS, STEP TO 0.1 KG

Réf. EKMWS01

KERN



CE
0297



SCALE RANGE 300 KGS, STEP TO 0.1 KG



Réf. **EKMWS01**

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0297	Barcode	4045761096684
Guarantee (month):	36	Customs Code	84231090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	80,0	Width packaging (cm):	95,0
Length (cm) :	115,5	Length packaging (cm) :	131,0
Height (cm) :	6,6	Height packaging (cm) :	12,0
Weight (kg) :	28,600	Weight packaging (kg) :	34,800
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:06 - This article is available on our website : www.gsh-med.fr with the following reference EKMWS01

All characteristics are given for information only and may be changed