

DESCHAMPS LIGATURE CARRIER RIGHT BLUNT, 21CM

Réf. IA05401



DESCHAMPS LIGATURE CARRIER RIGHT BLUNT, 21CM



Réf. IA05401

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040006260
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,8	Width packaging (cm):	
Length (cm) :	21,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,050	Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:42 - This article is available on our website : www.gsh-med.fr with the following reference IA05401

All characteristics are given for information only and may be changed