

STEVENS SCISSORS STRAIGHT S/S 10 CM

Réf. IC02410



STEVENS SCISSORS STRAIGHT S/S 10 CM



Réf. IC02410

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040008769
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,3	Width packaging (cm):	4,3
Length (cm) :	10,0	Length packaging (cm) :	10,0
Height (cm) :	0,3	Height packaging (cm) :	0,3
Weight (kg) :	0,019	Weight packaging (kg) :	0,020
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:19 - This article is available on our website : www.gsh-med.fr with the following reference IC02410

All characteristics are given for information only and may be changed