

ADSON HOOK SHARP 20 CM

Réf. IC19201



ADSON HOOK SHARP 20 CM



Réf. IC19201

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040010588
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,6	Width packaging (cm):	
Length (cm) :	20,2	Length packaging (cm) :	
Height (cm) :	0,6	Height packaging (cm) :	
Weight (kg) :	0,027	Weight packaging (kg) :	0,027
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference IC19201

All characteristics are given for information only and may be changed