

JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC

CA-MI

Réf. BCASPI77

JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC

CA-MI

Réf. BCASPI77

FEATURES			
Medical device	No		
Marquage	Barcode	3661040151076	
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,4	Width packaging (cm):	
Length (cm) :	27,0	Length packaging (cm) :	
Height (cm) :	9,3	Height packaging (cm) :	
Weight (kg) :	0,100	Weight packaging (kg) :	0,150
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:01 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI77

All characteristics are given for information only and may be changed