

KUNSCHER NAIL EXTRACTOR (WITHOUT HOOKS)

Réf. IC01416



KUNSCHER NAIL EXTRACTOR (WITHOUT HOOKS)



Réf. IC01416

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040008707
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	6,0
Length (cm) :		Length packaging (cm) :	20,5
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	0,350
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/06/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference IC01416

All characteristics are given for information only and may be changed