

**ADAPTOREU 100-240V~/50HZ/28MA AB: 01/12**  
**OUTPUT: 12V=/300MA/3,6VA**

Réf. EKB02



**ADAPTOREU 100-240V~/50HZ/28MA AB: 01/12**  
**OUTPUT: 12V=/300MA/3,6VA**

**Réf. EKB02**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040201009
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:16 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKB02

All characteristics are given for information only and may be changed