

CARDIOMATE 6 - ELECTROCARDIOGRAPH 6 CHANNELS

Réf. 260006



CARDIOMATE 6 - ELECTROCARDIOGRAPH 6 CHANNELS



Réf. 260006

ECG Cardiomate 3® and 6® are robust and reliable ECG. Featuring a 7" color screen tilt they allow viewing 12 lead simultaneously. Simple to use, recording and printing ECG are triggered with only one button. The Cardiomate 3® & 6® are Cardiomate equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer. The Cardiomate 3® & Cardiomate 6® are equipped with automatic interpretation and pacemaker detection. The internal battery life of 3 hours, their lightness and carrying handle make Cardiomate 3® & Cardiomate 6® the perfect companion in the office or traveling.

- Carrying handheld.
- Tilting screen : 7-inch color LCD display (resolution: 800x480 px)
- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM
- Autonomy :The operating life is 1080 ECG.
- Automatic interpretation and intervals calculation.
- Pacemaker detection.
- Network connection, USB & printer.
- Printing and rapid analysis.
- Acquisitions and display 12 simultaneous tracks.
- Alphanumeric keypad.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446025489
Guarantee (month):	24	Customs Code	9018110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	24,4	Width packaging (cm):	
Length (cm) :	31,0	Length packaging (cm) :	
Height (cm) :	6,5	Height packaging (cm) :	
Weight (kg) :	2,400	Weight packaging (kg) :	4,860
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	50
Humidity min	25	Humidity max	85



08/12/2022 12:38 - This article is available on our website : www.gsh-med.fr with the following reference 260006

All characteristics are given for information only and may be changed