

NON CHILL RING (BLACK) FOR CARDIOPRESTIGE II STETHOSCOPE SPENGLER

Réf. 507972



NON CHILL RING (BLACK) FOR CARDIOPRESTIGE II STETHOSCOPE SPENGLER



Réf. 507972

All characteristics are given for information only and may be changed

FEATURES			
Medical device	No		
Marquage		Barcode	3700446030315
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,001
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:48 - This article is available on our website : www.gsh-med.fr with the following reference 507972