

## SPO2 PEDIATRIC SILICON SENSOR FOR 4VS MONITOR



Réf. 551635



# SPO2 PEDIATRIC SILICON SENSOR FOR 4VS MONITOR



Réf. **551635**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446031220
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,083
SPECIFIC STORAGE CONDITION			
Temperature min	-40	Temperature max	55
Humidity min	5	Humidity max	95



31/03/2022 01:53 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 551635

All characteristics are given for information only and may be changed