

Réf. DSU231009



Réf. DSU231009

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	0120840935104046
Guarantee (month):	24		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	23,0
Length (cm) :			Length packaging (cm) :	25,0
Height (cm) :			Height packaging (cm) :	10,0
Weight (kg) :			Weight packaging (kg) :	0,766
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 01:58 - This article is available on our website : www.gsh-med.fr with the following reference DSU231009

All characteristics are given for information only and may be changed