

PEAK FLOWMETER

Réf. L36147



PEAK FLOWMETER



Réf. L36147

Un débitmètre offert pour 5 boîtes d'embouts

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0197	Barcode	3700446001780
Guarantee (month):	24	Customs Code	90261081
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,075
SPECIFIC STORAGE CONDITION			
Temperature min	1	Temperature max	50
Humidity min	10	Humidity max	95



31/03/2022 01:55 - This article is available on our website : www.gsh-med.fr with the following reference L36147

All characteristics are given for information only and may be changed