

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER

Réf. 221051



OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER



Réf. 221051

For SpO2 spot checks.
 Easy to use. Simply press the button to switch on.
 Reading comfort measures thanks to the rotatable TFT 2 colors screen :
 - Simultaneous display in less than 10 sec. :
 • Pulsed oxygen saturation (% SpO2),
 • Heart rate,
 • The pulse strength (bar graph),
 • Plethysmography curve.- Low battery indicator.
 - Customizable Display: 6 display modes.
 - 5 brightness levels.Switches off automatically after 5 seconds when not in use.
 Comes with wrist-strap, cover and 2 batteries.

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0123	Barcode	3700446025465
Guarantee (month):	24	Customs Code	90181910
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	7,3
Length (cm) :		Length packaging (cm) :	11,7
Height (cm) :		Height packaging (cm) :	4,9
Weight (kg) :		Weight packaging (kg) :	0,130
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/05/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference 221051

All characteristics are given for information only and may be changed