

## SUNTECH CUFF FOR SINGLE PATIENT NEONAT N°2 (4-8CM)

 SunTech Medical®

Réf. DSU231101



# SUNTECH CUFF FOR SINGLE PATIENT NEONAT N°2 (4-8CM)

Réf. DSU231101

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE		Barcode	0120840935103421
Guarantee (month):	0		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	13,0
Length (cm) :			Length packaging (cm) :	11,0
Height (cm) :			Height packaging (cm) :	16,0
Weight (kg) :			Weight packaging (kg) :	0,250
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 02:01 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DSU231101