

## BLACK BEZEL FOR SPENGLER ANEROID SPHYGMOMANOMETERS LIAN NANO AND VAQUEZ-LAUBRY NANO

Réf. 310163



# BLACK BEZEL FOR SPENGLER ANEROID SPHYGMOMANOMETERS LIAN NANO AND VAQUEZ-LAUBRY NANO



**Réf. 310163**

Spare part SPENGLER, Black bezel for SPENGLER aneroid sphygmomanometers Lian Nano and Vaquez-Laubry Nano

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3700446028626
Guarantee (month):	60	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:42 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 310163

All characteristics are given for information only and may be changed