

THERMOMETRES TEMPO 10 FLEX GREY

Réf. SP1648



THERMOMETRES TEMPO 10 FLEX GREY



Réf. **SP1648**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446041618
Guarantee (month):	24	Customs Code	9025192090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	6,0
Length (cm) :		Length packaging (cm) :	15,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	0,052
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:57 - This article is available on our website : www.gsh-med.fr with the following reference SP1648

All characteristics are given for information only and may be changed