

OPD CUFF ADULT 23-33CM METAL MALE BAYONET (BOX OF 5)

Réf. DSU231113

SunTech Medical®



OPD CUFF ADULT 23-33CM METAL MALE BAYONET (BOX OF 5)

Réf. DSU231113

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661041008812
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	9,0
Height (cm) :		Height packaging (cm) :	18,0
Weight (kg) :		Weight packaging (kg) :	0,450
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:56 - This article is available on our website : www.gsh-med.fr with the following reference DSU231113