

## BABY TO TODDLER ELECTRONIC SCALE + BAG

Réf. EHANDI3S



## BABY TO TODDLER ELECTRONIC SCALE + BAG



Réf. EHANDI3S

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3361040000021
Guarantee (month):	12	Customs Code	84231090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	34,8
Length (cm) :		Length packaging (cm) :	63,5
Height (cm) :		Height packaging (cm) :	16,9
Weight (kg) :		Weight packaging (kg) :	1,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:12 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EHANDI3S

All characteristics are given for information only and may be changed

## BABY TO TODDLER ELECTRONIC SCALE + BAG



Réf. EHANDI3S

### COMPOSITION

Code	Descriptive	Quantity
EHANDI3	BABY TO TODDLER ELECTRONIC SCALE	1
EHSACTA3	TRANSPORT BAG FOR ELECTRONIC BABY WEIGHT SCALE HOLTEX not Sold alone	1