

Réf. DSU230700



Réf. DSU230700

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0413	Barcode	10840935105701
Guarantee (month):	24	Customs Code	9018.90.5040
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	17,4	Width packaging (cm):	
Length (cm) :	24,0	Length packaging (cm) :	
Height (cm) :	11,5	Height packaging (cm) :	
Weight (kg) :	1,680	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:41 - This article is available on our website : www.gsh-med.fr with the following reference DSU230700