

ALUMINIUM ALLOY STAIR STRETCHER -W 46 CM - 4 WHEELS, ORANGE

Réf. SBRAN15_O



ALUMINIUM ALLOY STAIR STRETCHER -W 46 CM - 4 WHEELS, ORANGE



Réf. **SBRAN15_0**

- High strength aluminium alloy structure
- Equipped with armrests
- Ideal for transferring a patient, climbing up and down the stairs
- Back wheels are fixed and have brakes, the front two are swivel wheels, no brake. This allowing movement on most surfaces with a minimum of effort.
- The safety of patients during the transport process is ensured by 2 restraining straps with quick release buckle
- At the front of the stretcher, two retractable arms
- Dimensions (L x W x H):
90 X 51 x 91 cm, unfolded;
90 x 17 x 59 cm, folded
- Weight: 10.5 kg
- Maximum load carried: 159 kg
- Rigid seat
- Single colour: orange

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040094694
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	60,0	Width packaging (cm):	64,0
Length (cm) :	84,0	Length packaging (cm) :	10,0
Height (cm) :	98,0	Height packaging (cm) :	20,0
Weight (kg) :	14,000	Weight packaging (kg) :	16,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:21 - This article is available on our website : www.gsh-med.fr with the following reference SBRAN15_0