

MALE CONNECTOR (FOR VACUUM PORT) FOR ASKIR 20 PACK OF 5'S

Réf. BCASPI164

CA-MI

GARANTIE
1
AN
TVA
GARANTIE

CE



MALE CONNECTOR (FOR VACUUM PORT) FOR ASKIR 20 PACK OF 5'S

CA-MI

Réf. BCASPI164

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3661040160375
Guarantee (month):	12	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,010	Weight packaging (kg) :	0,010
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:08 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI164

All characteristics are given for information only and may be changed