

ROLL OF SILICONE TUBE Ø8X14MMX10M ASK 36BR/118/C30/C30BR/NHOSPIVAC

Réf. BCASPI145

CA-MI



CE
0123



ROLL OF SILICONE TUBE Ø8X14MMX10M ASK 36BR/118/C30/C30BR/NHOSPIVAC



Réf. BCASPI145

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | No | | |
| Marquage | CE0123 | Barcode | 3661040163925 |
| Guarantee (month): | 12 | Customs Code | |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | 39,0 |
| Length (cm) : | | Length packaging (cm) : | 39,0 |
| Height (cm) : | | Height packaging (cm) : | 8,0 |
| Weight (kg) : | | Weight packaging (kg) : | 1,422 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | -25 | Temperature max | 70 |
| Humidity min | 1 | Humidity max | 85 |



30/03/2022 23:36 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI145