

PEARL MAGENTA STAINLESS STEEL SINGLE TUBE BINAURALS

Réf. 507898



PEARL MAGENTA STAINLESS STEEL SINGLE TUBE BINAURALS

Réf. **507898**

Spare part for Magister® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446003920
Guarantee (month):	60		Customs Code	90181990
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,080	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:46 - This article is available on our website : www.gsh-med.fr with the following reference 507898

All characteristics are given for information only and may be changed