

## BLACK CHROME PLATED DUAL TUBE BINAURALS

Réf. 507920



# BLACK CHROME PLATED DUAL TUBE BINAURALS

Réf. **507920**

Spare part for Cardio Sonic® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446009588
Guarantee (month):	0		Customs Code	90181990
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,095
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 01:47 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507920

All characteristics are given for information only and may be changed