

NON CHILL RING (GREEN) FOR PULSE II PEADIATRIC STETHOSCOPE SPENGLER



Réf. 507983

NON CHILL RING (GREEN) FOR PULSE II PEADIATRIC STETHOSCOPE SPENGLER



Réf. **507983**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446030421
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,001
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:49 - This article is available on our website : www.gsh-med.fr with the following reference 507983

All characteristics are given for information only and may be changed